

# 高效率金屬移除劑 Metal Scavenger



Uni-Onward  
友和電子報

產品介紹

產品選擇指南

操作示範影片

## SiliaMetS® Metal Scavengers

SiliCycle, THE world leader  
for more than 12 years  
in metal scavenging solutions!



### SiliaMetS® 金屬移除劑的特性&優點

1. 操作方便
2. 不污染
3. 可移除多種金屬催化劑
4. 室溫即可使用
5. 適用各種水溶液&有機溶劑PH2-12



加入適量的移除劑在溶液中攪拌後簡單過濾即可得到清澈solution

Product item	Silia MetS	Ag	Ca	Cd	Co	Cr	Cs	Cu	Fe	Hg	Ir	La	Li	Mg	Ni	Os	Pb	Pd (II)	Pd (0)	Pt	Rh (1)	Rh (II)	Rh (III)	Ru (II)	Sc	Sn	W	Zn	
SY-R51030B	<u>Thiol</u>	★						☆		★	☆						★	☆	★	★		☆	☆	☆	★	☆	☆		
SY-R69530B	<u>Thiourea</u>	☆						☆	☆								☆		★	★		☆	☆	☆		☆	☆		
SY-R80530B	<u>Cysteine</u>		☆	★		☆	☆	☆	★		☆			☆		★		☆	☆	☆	☆	☆		★	★	★		☆	
SY-R79030B	<u>DMT</u>			☆	☆			☆	☆		★				★	★		★	★	★	★	★	★	★	☆			☆	
SY-R52030B	<u>Amine</u>			★	☆	★		☆	☆	☆							☆	☆		★	★	★					☆	☆	
SY-R49030B	<u>Diamine</u>			☆	☆	★		☆	☆	☆					☆		☆	★	★	★				☆			☆	☆	
SY-R48030B	<u>Triamine</u>				☆	★		☆	☆						☆		★	★	★	★				☆			☆	☆	
SY-R79230B	<u>Imidazole</u>			★	★	☆		★	★						★	★		☆	☆		☆	☆					★	★	
SY-R69030B	<u>TAAcOH</u>				★	☆	☆		☆						★	★		☆	☆		☆	☆			★	☆			
SY-R69230B	<u>TAAcONa</u>		★	★		☆	★	★	★		★	★	★	★	★	★		☆	☆		☆	☆	★	☆	★	★		☆	

★ Preferred Scavengers / ☆ Scavenges



1. 產品介紹:如附件DM

2. 產品選擇指南連結網址

SiliaMetS® Metal Scavenger Selection Table

<http://www.silicycle.com/products/metal-scavengers/siliabond-metal-scavengers-selection-guide>

3. 操作示範影片

<https://www.youtube.com/watch?v=vsTHXboxA4o>

# 為您的食品安全把關! Mycotoxin Analysis 黴菌毒素分析



黴菌毒素是真菌的二次代謝產物，其生產之二次代謝物含有對人、畜有害的成分。黴菌毒素其英文名中之“myco”是希臘語 mykes (黴菌) 之意，“toxin”則為拉丁語 toxicum (毒) 之意。食物和飼料在生長、採收、運送以及儲藏過程都可能遭遇黴菌毒素污染，黴菌毒素的存在可能會影響人類和動物的健康、最常見的黴菌毒素污染源：如黃麴毒素 (aflatoxin B1, B2, G1 and G2, aflatoxin M1) · 赭麴毒素A (OTA) · 單端孢霉烯族毒素類 (Deoxynivalenol –DON, nivalenol, T-2 toxin, HT-2 toxin) 和伏馬毒素B1和B2等毒素。

SIMGA 和 MN 兩大知名廠牌提供您 MYCOTOXIN ANALYSIS WORKFLOW 最完整的相關產品，一起為您的食品安全把關。

## MYCOTOXIN ANALYSIS WORKFLOW & PRODUCTS



### Supel Tox Products SUPELCO

SPE Cartridge	Analyte(s)	Qty.	Cat. No.
Supel Tox AflaZea	Aflatoxin B <sub>1</sub> , B <sub>2</sub> , G <sub>1</sub> and G <sub>2</sub> , and Zearalenone	30	55314-U
Supel Tox DON	Deoxynivalenol (DON)	30	55316-U
Supel Tox Tricho	Trichothecenes (Type A and B)	30	55308-U
Supel Tox TrichoBind	Trichothecenes (Type A and B)	25	55307-U
Supel Tox FumoniBind	Fumonisin (B <sub>1</sub> and B <sub>2</sub> )	25	55315-U
Supel Tox OchraBind	Ochratoxin A	25	55318-U

[更多產品訊息](#)

<http://www.sigmaaldrich.com/catalog/search?term=Supel%E2%84%A2+Tox+SPE+Cartridges&interface=All&N=0&mode=match%20partialmax&lang=en&region=TW&focus=product>


### SPE Products MACHEREY-NAGEL

Step	Product	Pack of	REF
SPE	CHROMABOND® Carbon/Alox/Celite®, 6 mL, 500 mg (on request)	30	730584
	CHROMABOND® vacuum manifold for 12 columns	1	730150

[更多產品訊息](#)

<http://www.mn-net.com/tabid/12037/language/en-US/default.aspx>

Cleanup of Aflatoxins in Peanut Paste: Immunoaffinity Versus Supel™ Tox (n = 3)



Supel Tox SPE 新型黴菌毒素專用高速固相萃取管相對於傳統萃取親合管，具有10倍之萃取速度及更高的再現性。

### Supel Tox SPE Cartridges Mycotoxins

Supel™ Tox SPE Cartridges  
Fast and Simple Cleanup for Mycotoxin Analysis

Think about food  
Care about safety

Vol. 6  
Mycotoxins

MACHEREY-NAGEL

Matrix Processing

Isolation and Concentration

Derivatization

Chromatographic Separation

Quantification

## Mycotoxin Standard

**Fluka**  
Analytical

- Mycotoxin reference materials (Neats)
- Mycotoxin single component standard solutions
- Mycotoxin mixture standard solutions for multi-analyte detection
- Dried down mycotoxin standards (RMs)
- Certified matrix reference materials for mycotoxin analysis (CRMs)
- Mycotoxin solutions (CRMs) from IRMM

更多產品訊息請點選電子期刊

## Mycotoxin Standard



## Mycotoxin Analysis

